

Food and Environmental Proficiency Testing Unit

The final date for return of results is indicated on the accompanying request/report form

Contact details:

The Organisers - FEPTU UK Health Security Agency 61 Colindale Avenue, London, NW9 5EQ, UK. Tel: +44 (0) 20 8327 7119

e-mail: foodeqa@ukhsa.gov.uk



Results should be returned on-line using the web-based reporting system. Go to: www.phe-eqa.org.uk. You will need your log-in details for this process. Please contact us if you experience problems with the on-line system.

0

European Food Microbiology Legislation - Food and Environmental Proficiency Testing Unit

European Food Microbiology Legislation Scheme – Instruction Sheet

Each distribution consists of **three** LENTICULE® discs in screw-cap plastic vials (with desiccant). The LENTICULE discs require reconstitution by a process of re-hydration and dispersion prior to examination.

Refer also to the Safety Data Sheet: Safety data sheet - LENTICULE

For a video on sample processing: Processing samples - European Food Law

Storage:

- Store the samples at -20 ± 5°C on receipt.
- Allow the LENTICULE discs to reach ambient temperature (5 10 minutes) before re-hydrating in diluent.

Reconstitution:

- Open the sample container and transfer the LENTICULE disc into 100mL buffered peptone water (BPW).
- Leave at ambient temperature for a minimum of 10 minutes to hydrate.
- Ensure that the LENTICULE disc has completely hydrated and disperse the inoculum by inverting approximately 30 times.

Examination:

- Each reconstituted sample is equivalent to a 100g food sample.
- Undertake sample examination between 30 45 minutes after reconstitution.
- Examine in accordance with routine procedures.

Return of Results:

- Return your results on-line using the web-based reporting system.
- Return results to the organisers before the deadline shown on the request/report form.
- Results returned after that date will not be assessed.

FEPTU382.06 Page 1 of 1