

Food and Environmental Proficiency Testing Unit

The final date for return of results is indicated on the accompanying request/report form

Contact details:

The Organisers - FEPTU UK Health Security Agency 61 Colindale Avenue, London, NW9 5EQ, UK.

Tel: +44 (0) 20 8327 7119 e-mail: foodeqa@ukhsa.gov.uk

Results to be returned via an online survey as per email sent to your laboratory

Norovirus and Hepatitis A - Food and Environmental Proficiency Testing Unit

Norovirus and Hepatitis A Virus Scheme - Instruction Sheet

This scheme consists of **two** LENTICULE® discs in screw-cap plastic vials (with desiccant). The LENTICULE discs require reconstitution by a process of re-hydration and dispersion prior to examination.

Refer also to the Safety Data Sheet: Safety data sheet - LENTICULE

Storage:

- Store the samples at -20 ± 5°C on receipt.
- Allow the LENTICULE discs to reach ambient temperature (5 10 minutes) before rehydrating in diluent.

Reconstitution:

- Open the plastic tube and invert it to transfer the LENTICULE disc into a new clean tube.
- Add 975 µL of sterile phosphate buffered saline (PBS) to the LENTICULE disc.
- Leave at ambient temperature for a minimum of 10 minutes to rehydrate. Ensure that the LENTICULE disc has dissolved completely before proceeding, this becomes the sample.
- Disperse the inoculum by inverting approximately 30 times.
- Examine samples between 30 45 minutes after reconstitution.

It is not necessary to include a specific virus extraction step – the samples can be subjected directly to RNA extraction. Analyse each sample for all target viruses. If using a quantitative real-time method, quantify each test sample.

Examination:

- Examine sample using your routine molecular methods.
- Detection and quantification (if routinely done) of viruses is required for both samples.

Return of Results:

- Return your results via an online survey as per email sent to your laboratory.
- Return results to the organisers before the deadline shown on the request/report form.

FEPTU830.05 Page 1 of 1