

# Food and Environmental Proficiency Testing Unit

Laboratory identification no. (check): <Lab No>  
 Dispatch date: 03 October 2022  
 Final date for return of results: 28 October 2022

**Contact details:**

The Organisers - FEPTU  
 UK Health Security Agency  
 61 Colindale Avenue  
 London, NW9 5EQ, UK.  
 Fax: +44 (0) 20 8200 8264  
 Tel: +44 (0) 20 8327 7119  
 e-mail: foodeqa@ukhsa.gov.uk



[www.gov.uk/government/collections/external-quality-assessment-ega-and-proficiency-testing-pt-for-food-water-and-environmental-microbiology#dialysis-water-scheme](http://www.gov.uk/government/collections/external-quality-assessment-ega-and-proficiency-testing-pt-for-food-water-and-environmental-microbiology#dialysis-water-scheme)

## Dialysis Water Scheme Request/Report Form

<b>Distribution No: DW28</b>	<b>Sample numbers: DW28A and DW28B</b>
Download the sample instruction sheet. A video on sample processing is located further down the webpage:	<a href="http://www.gov.uk/government/publications/dialysis-water-scheme-sample-instruction-sheet">www.gov.uk/government/publications/dialysis-water-scheme-sample-instruction-sheet</a>
Download the safety data sheet:	<a href="http://www.gov.uk/government/publications/safety-data-sheet-for-lenticules">www.gov.uk/government/publications/safety-data-sheet-for-lenticules</a>
<b><i>If you cannot examine any of these return your results as 'Not examined'</i></b>	
Sample type:	Water sample collected three hours ago from a sample port in the distal segment of the loop of the dialysis machine.

**Results:**

Sample Number	Enumeration	Not examined	Analyst	Conclusion reported to the sender (Please tick any that apply)
DW28A	TVC 17°C - 23°C for 7 days per mL	<input type="checkbox"/>	Nominated result:	Satisfactory/Acceptable <input type="checkbox"/>
				Unsatisfactory/Unacceptable <input type="checkbox"/>
				Borderline <input type="checkbox"/>
			Result 2:	Breaches action level <input type="checkbox"/>
				Other, please state <input type="checkbox"/>
				Meets the maximum admitted value <input type="checkbox"/>
			Result 3:	Satisfactory/Triggers action level <input type="checkbox"/>
				Not routinely reported <input type="checkbox"/>

# Food and Environmental Proficiency Testing Unit

Laboratory identification no. (check): <Lab No>

Sample Number	Enumeration	Not examined	Analyst	Conclusion reported to the sender (Please tick any that apply)	
DW28B	TVC 17°C - 23°C for 7 days per mL	<input type="checkbox"/>	Nominated result:	Satisfactory/Acceptable	<input type="checkbox"/>
				Unsatisfactory/Unacceptable	<input type="checkbox"/>
				Borderline	<input type="checkbox"/>
			Result 2:	Breaches action level	<input type="checkbox"/>
				Other, please state	<input type="checkbox"/>
				Meets the maximum admitted value	<input type="checkbox"/>
			Result 3:	Satisfactory/Triggers action level	<input type="checkbox"/>
				Not routinely reported	<input type="checkbox"/>

Please indicate the method used	<input type="checkbox"/>	Pour
	<input type="checkbox"/>	Membrane filtration
	<input type="checkbox"/>	Spread
	<input type="checkbox"/>	IDEXX EasyDisc YEA
	<input type="checkbox"/>	Not stated
	<input type="checkbox"/>	Other
Please indicate the media used	<input type="checkbox"/>	Tryptone glucose extract agar (TGEA)
	<input type="checkbox"/>	Yeast extract agar (YEA)
	<input type="checkbox"/>	Reasoner's 2A agar (R2A)
	<input type="checkbox"/>	Tryptone soya agar (TSA)
	<input type="checkbox"/>	Other
Standard/Guideline followed	<input type="checkbox"/>	ISO 23500-3:2019
	<input type="checkbox"/>	ISO 23500-5:2019
	<input type="checkbox"/>	ISO 6222:1999
	<input type="checkbox"/>	European Best Practice Guidelines for Haemodialysis
	<input type="checkbox"/>	UK Renal Association guideline
	<input type="checkbox"/>	Public Health England guideline
	<input type="checkbox"/>	European Pharmacopoeia 5.0
Additional comments	<input type="checkbox"/>	Other, please state

Authorised by:

Date reported: